

Product Datasheet

Recombinant Mouse Interleukin-2/IL-2 (orb1906494)

Description	Recombinant Mouse Interleukin-2 is produced by our E.coli expression system and the target gene enco...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Sodium Citrate, 0.2%Tween 80,pH3.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Ala21-Gln169
Source	E.coli
Uniprot ID	P04351
Expiration Date	6 months from date of receipt.