

Product Datasheet

Recombinant Mouse IL-1 Receptor Antagonist Protein/IL-1ra/IL-1RA (orb1906495)

Description	Recombinant Mouse Interleukin-1 Receptor Antagonist Protein is produced by our E.coli expression sys...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Arg27-Gln178
Source	E.coli
Uniprot ID	P25085
Expiration Date	6 months from date of receipt.