

Product Datasheet

Recombinant Mouse HVEM/TNFRSF14/TR2/CD270 (C-Fc) (orb1906496)

Description	Recombinant Mouse Herpesvirus Entry Mediator is produced by our Mammalian expression system and the ...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Gln39-Val207
Source	Human Cells
Uniprot ID	Q80WM9
Expiration Date	6 months from date of receipt.