

Product Datasheet

Recombinant Human IL-2 Receptor Subunit α /IL-2RA/CD25 (C-Fc) (orb1906498)

Description	Recombinant Human Interleukin-2 Receptor subunit alpha is produced by our Mammalian expression syste...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Protein Sequence	Glu22-Cys213
Source	Human Cells
Uniprot ID	P01589
Expiration Date	6 months from date of receipt.