

Product Datasheet

Recombinant Mouse Apolipoprotein E/ApoE (C-6His) (orb1906500)

Description	Recombinant Mouse Apolipoprotein E is produced by our Mammalian expression system and the target gen...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Protein Sequence	Glu19-Gln311
Source	Human Cells
Uniprot ID	P08226
Expiration Date	6 months from date of receipt.