

Product Datasheet

Recombinant Human Actin-related protein 2/3 complex subunit 3/ARPC3 (orb1906631)

Description	Recombinant Human Actin-related protein 2/3 complex subunit 3/ARPC3 is produced by E. coli. The target...
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-FLAG
Note	For research use only
Protein Sequence	P2-Q178
Source	E. coli
Uniprot ID	O15145
Expiration Date	6 months from date of receipt.