

Product Datasheet

Recombinant Human L-lactate dehydrogenase A-like 6A/LDHAL6A (orb1906687)

Description	Recombinant Human L-lactate dehydrogenase A-like 6A/LDHAL6A is produced by E. coli. The target gene e...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	P210-L332
Source	E. coli
Uniprot ID	Q6ZMR3
Expiration Date	6 months from date of receipt.