

Product Datasheet

Recombinant Human Prostaglandin reductase 2/PTGR2 (orb1906741)

Description	Recombinant Human Prostaglandin reductase 2/PTGR2 is produced by E. coli. The target gene encoding l2...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-FLAG
Note	For research use only
Protein Sequence	I2-L351
Source	E. coli
Uniprot ID	Q8N8N7
Expiration Date	6 months from date of receipt.