

Product Datasheet

Recombinant Human UCHL5 Protein (orb1906742)

Description	Recombinant Human UCHL5 Protein is produced by E. coli. The target gene encoding T2-K315 is expressed ...
Form/Appearance	Supplied as a 0.2 µm filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-FLAG
Note	For research use only
Protein Sequence	T2-K315
Source	E. coli
Uniprot ID	Q9Y5K5
Expiration Date	6 months from date of receipt.