

Product Datasheet

Recombinant Human VPS29 Protein (orb1906746)

Description	Recombinant Human VPS29 Protein is produced by E. coli. The target gene encoding N68-P182 is expressed...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-FLAG
Note	For research use only
Protein Sequence	N68-P182
Source	E. coli
Uniprot ID	Q9UBQ0
Expiration Date	6 months from date of receipt.