

Product Datasheet

Recombinant Human Mitotic checkpoint protein BUB3 (orb1906792)

Description	Recombinant Human Mitotic checkpoint protein BUB3 is produced by E. coli. The target gene encoding T2-...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 10% Glycerol, 300mM NaCl, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	T2-Y141
Source	E. coli
Uniprot ID	O43684
Expiration Date	6 months from date of receipt.