



Product Datasheet

CREB3L1 antibody (orb19068)


Description	Goat polyclonal antibody to CREB3L1
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	OASIS / CREB3L1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in lysates of cell line A549 (calculated MW of 57.0kDa according to NP_443086.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DFSHSKEWFHDRD
Clonality	Polyclonal
MW	57
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 lysate using...



Immunohistochemical staining of Human Sk...



Immunohistochemical staining of Human Li...