

Product Datasheet

Recombinant Human Small G protein signaling modulator 3/SGSM3 (orb1906913)

Description	Recombinant Human Small G protein signaling modulator 3/SGSM3 is produced by E. coli. The target gen...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	F529-E724
Source	E. coli
Uniprot ID	Q96HU1
Expiration Date	6 months from date of receipt.