

Product Datasheet

Recombinant Human MOB-like protein phocein/MOB4 (orb1906916)

Description	Recombinant Human MOB-like protein phocein/MOB4 is produced by E. coli. The target gene encoding V2-A2...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	V2-A225
Source	E. coli
Uniprot ID	Q9Y3A3
Expiration Date	6 months from date of receipt.