

Product Datasheet

Recombinant Human Ras association domain-containing protein 3/RASSF3 (orb1906924)

Description	Recombinant Human Ras association domain-containing protein 3/RASSF3 is produced by E. coli. The targ...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	K84-D238
Source	E. coli
Uniprot ID	Q86WH2
Expiration Date	6 months from date of receipt.