

Product Datasheet

Recombinant Human Retinoid-Binding Protein 7 (orb1906950)

Description	Recombinant Human Retinoid-binding Protein 7 is produced by our E.coli expression system and the tar...
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Ala134
Source	E. coli
Uniprot ID	Q96R05
Expiration Date	6 months from date of receipt.