

Product Datasheet

Recombinant Human SSU72 Protein (orb1906954)

Description	Recombinant Human SSU72 Protein is produced by E. coli. The target gene encoding P2-Y194 is expressed ...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10% Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	P2-Y194
Source	E. coli
Uniprot ID	Q9NP77
Expiration Date	6 months from date of receipt.