

Product Datasheet

Recombinant Human E3 ubiquitin-protein ligase RNF114 (orb1907029)

Description	MADE TO ORDER Recombinant Human E3 ubiquitin-protein ligase RNF114 is produced by E.coli expressio...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	E119-Q228
Source	E. coli
Uniprot ID	Q9Y508
Expiration Date	6 months from date of receipt.