

Product Datasheet

Recombinant Human PSMD13 Protein (orb1907064)

Description	MADE TO ORDER Recombinant Human PSMD13 Protein is produced by E.coli expression system and the tar...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	K2-W232
Source	E. coli
Uniprot ID	Q9UNM6
Expiration Date	6 months from date of receipt.