

Product Datasheet

Recombinant Human Glycogen phosphorylase, liver form (orb1907129)

Description	MADE TO ORDER Recombinant Human Glycogen phosphorylase, liver form is produced by E.coli expressio...
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, 300mM NaCl, 10%Glycerol, pH 7.4
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	V637-N847
Source	E. coli
Uniprot ID	P06737
Expiration Date	6 months from date of receipt.