

Product Datasheet

Recombinant Human D-Tyrosyl-tRNA(Tyr) Deacylase 1/DTD1 (C-6His) (orb1907181)

Description	Recombinant Human D-Tyrosyl-tRNA(Tyr) Deacylase 1 is produced by our mammalian expression system and...
Form/Appearance	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	Met1-Pro209
Source	Human Cells
Uniprot ID	Q8TEA8
Expiration Date	6 months from date of receipt.