

## Product Datasheet

### KAAG1 Antibody (orb1907256)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | KAAG1 Antibody                                                                   |
| <b>Reactivity</b>          | Human                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                     |
| <b>Tested Applications</b> | ELISA                                                                            |
| <b>Immunogen</b>           | Recombinant Human KAAG1 protein                                                  |
| <b>Target</b>              | KAAG1                                                                            |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| <b>Form/Appearance</b>     | Liquid                                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | hIgG1                                                                            |
| <b>Clonality</b>           | Monoclonal                                                                       |
| <b>Uniprot ID</b>          | <b>Q9UBP8</b>                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |