

Product Datasheet

ROBO1 Antibody (orb1907265)

Description	ROBO1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Immunogen	Recombinant Human ROBO1 protein
Target	ROBO1
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	Q9Y6N7
Expiration Date	12 months from date of receipt.