

Product Datasheet

IGF1R Antibody (orb1907283)

Description	IGF1R Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Recombinant Human IGF1R protein
Target	IGF1R
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P08069
Expiration Date	12 months from date of receipt.