

Product Datasheet

HIF1A Antibody (orb1907292)

Description	HIF1A Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Immunogen	Recombinant Human HIF1A protein
Target	HIF1A
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	Q16665
Expiration Date	12 months from date of receipt.