


Product Datasheet

VGLL4 antibody (orb1908)

Description	VGLL4 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human VGLL4 (201-290/290 aa)
Target	VGLL4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	Q14135
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



WB analysis
of Ovary
(Mouse)
Lysate at 4...