



Product Datasheet

MEFV antibody (orb19088)

Description	Goat polyclonal antibody to MEFV
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	MEFV
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 80kDa band observed in Human Peripheral Blood Leukocyte lysate. Recommended concentration: 1-3 µg/ml during overnight incubation.
Protein Sequence	EHLKKLRKSGEEQ
Clonality	Polyclonal
MW	86.4; 50.9
Entrez	4210
Dilution Range	ELISA: 1:8000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A431 cell...



Flow Cytometry analysis of A431 cells of...