

Product Datasheet

GOLGA6B antibody (orb19090)

Description

Goat polyclonal antibody to GOLGA6B

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

GOLGA6 / GLP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:1000. WB: Preliminary experiments gave an approx 45kDa band in human liver lysates after 1 µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 79.9kDa according to NP_001033729.1. The 45kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

DLPKEKADGTEQ

Clonality

Polyclonal

MW


79.9

Entrez
[55889](#), [342096](#)
Dilution Range

ELISA: 1:1000, WB: 1 µg/ml, IHC-P: 3-6µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of Human Ce...