

Product Datasheet

Mouse Serine Protease Inhibitor A3N (Serpina3n) ELISA Kit (orb1909607)

Description	Mouse Serine Protease Inhibitor A3N (Serpina3n) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q91WP6
Sensitivity	0.124ng/mL
Expiration Date	Please enquire.