

Product Datasheet

Human Secretory Carrier-Associated Membrane Protein 3(SCAMP3) ELISA Kit (orb1909608)

Description	Human Secretory Carrier-Associated Membrane Protein 3(SCAMP3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O14828
Sensitivity	0.062ng/mL
Expiration Date	Please enquire.