

Product Datasheet

Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit (orb1909609)

Description	Rat Runt Related Transcription Factor 2 (RUNX2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9Z2J9
Sensitivity	0.071ng/mL
Expiration Date	Please enquire.