

Product Datasheet

Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit (orb1909622)

Description	Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9HC73
Sensitivity	0.052ng/mL
Expiration Date	Please enquire.