

Product Datasheet

Rat Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit (orb1909633)

Description	Rat Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q64666
Sensitivity	0.137ng/mL
Expiration Date	Please enquire.