

Product Datasheet

Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit (orb1909645)

Description	Mouse Low Density Lipoprotein Receptor Related Protein 6 (LRP6) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O88572
Sensitivity	0.067ng/mL
Expiration Date	Please enquire.