

## Product Datasheet

# Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit (orb1909653)

|                        |                                                                                          |
|------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>     | Human Cytokine Receptor Like Factor 2 (CRLF2) ELISA Kit                                  |
| <b>Reactivity</b>      | Human                                                                                    |
| <b>Range</b>           | 0.156-10ng/mL                                                                            |
| <b>Storage</b>         | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>            | For research use only                                                                    |
| <b>Sample Types</b>    | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Uniprot ID</b>      | <b>Q9HC73</b>                                                                            |
| <b>Sensitivity</b>     | 0.052ng/mL                                                                               |
| <b>Expiration Date</b> | Please enquire.                                                                          |