

Product Datasheet

Human Cyclic ADP Ribose Hydrolase (cADPRH) ELISA Kit (orb1909659)

Description	Human Cyclic ADP Ribose Hydrolase (cADPRH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P28907
Sensitivity	0.068ng/mL
Expiration Date	Please enquire.