

Product Datasheet

Human B-Cell CLL/Lymphoma 6(Bcl6) ELISA Kit (orb1909662)

Description	Human B-Cell CLL/Lymphoma 6(Bcl6) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P41182
Sensitivity	0.112ng/mL
Expiration Date	Please enquire.