

Product Datasheet

Rat Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit (orb1909664)

Description	Rat Arylalkylamine-N-Acetyltransferase (AANAT) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q64666
Sensitivity	0.137ng/mL
Expiration Date	Please enquire.