

Product Datasheet

Rat Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit (orb1909665)

Description	Rat Alpha-2-Plasmin Inhibitor (a2PI) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma and other biological fluids.
Sensitivity	0.115ng/mL
Expiration Date	Please enquire.