

Product Datasheet

Fish Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) ELISA Kit (orb1909698)

Description	Fish Tyrosine Kinase With Immunoglobulin Like And EGF Like Domains Protein 1 (Tie1) ELISA Kit
Reactivity	Fish
Range	7.81-500ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	073792
Sensitivity	3.2ng/mL
Expiration Date	Please enquire.