

Product Datasheet

Fish Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit (orb1909702)

Description	Fish Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit
Reactivity	Fish
Range	0.156-10ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.086ng/mL
Expiration Date	Please enquire.