

Product Datasheet

Fish Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit (orb1909703)

Description	Fish Insulin Like Growth Factor Binding Protein 1 (IGFBP1) ELISA Kit
Reactivity	Fish
Range	0.156-10ng/mL
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.086ng/mL
Expiration Date	Please enquire.