

Product Datasheet

Bovine Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit (orb1909704)

Description	Bovine Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Storage	All reagents should be stored at -20°C upon receipt
Note	For research use only
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00442
Sensitivity	0.120ng/mL
Expiration Date	Please enquire.