

Product Datasheet

**Human LTBR / TNFRSF3 Protein, His, (MALS
verified) (orb1909735)**

Description

Human LTBR Protein, His, is

(HEK293). It contains AA Gln 31 - Met 227...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LTBR

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 25.3 kDa. The protein migrates as 33-35 kDa and 37-40 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Human LTBR Protein, His, is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Met 227 (Accession # P36941-1).

MW

25.3 kDa

Uniprot ID
P36941
Expiration Date

6 months from date of receipt.