

Product Datasheet

Human CCL8 Protein, His Tag (orb1909738)

Description

Human CCL8 Protein, His Tag is expressed

	Pro 99 (Accessi...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CCL8
Preservatives	Powder
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 10.9 kDa. The protein migrates as 10-11 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).
Purity	90%
Source	Human CCL8 Protein, His Tag is expressed from E. coli cells. It contains AA Gln 24 - Pro 99 (Accession # P80075).
MW	10.9 kDa
Uniprot ID	P80075
Expiration Date	6 months from date of receipt.