



Product Datasheet

CD160 Antibody (FITC) (orb1909788)


Description	CD160 is a cell surface glycoprotein of immunoglobulin superfamily, which functions as a costimulato...
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Human NK cell line YT2C2
Target	CD160
Preservatives	Stabilizing Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgM kappa
Clonality	Monoclonal
Uniprot ID	O95971
Entrez	11126
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD160
positive
CD56 ...



Flow
cytometry
surface
staining
pattern ...