



Product Datasheet

C5aR2 Antibody (APC) (orb1909789)

Description	C5aR2, also known as C5L2, is one of two receptors for C5a (anaphylatoxin). It is coexpressed with C...
Reactivity	Human, Primate
Conjugation	APC
Tested Applications	FC
Immunogen	L1.2 cells transfected with human C5aR2
Target	C5aR2
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Uniprot ID	Q9P296
Entrez	27202
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocy...



Flow cytometry surface staining pattern ...