


Product Datasheet

CD105 Antibody (APC) (orb1909791)

Description	CD105 (endoglin) is a homodimeric transmembrane glycoprotein serving in presence of TGFbetaR-2 as a ...
Reactivity	Human, Porcine
Conjugation	APC
Tested Applications	FC
Immunogen	Recombinant Vaccinia virus containing the human CD105 (L-isoform) cDNA.
Target	CD105
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Uniprot ID	P17813
Entrez	2022
Expiration Date	12 months from date of receipt.



Separation of HUVEC cells stained using ...