



Product Datasheet

PCK2 antibody (orb19100)

Description	Goat polyclonal antibody to PCK2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	PCK2 / PEPCK-M
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 70kDa band observed in Human Liver lysates and in lysates of cell line HepG2 (calculated MW of 70.7kDa according to NP_004554.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KPWKPGDKEPCAH
Clonality	Polyclonal
MW	70.7; 56.0
Entrez	5106
Dilution Range	ELISA: 1:4000, WB: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysate US...



Western blot analysis of Human Liver lys...