

Product Datasheet

DAZL antibody (orb19110)

Description

Goat polyclonal antibody to DAZL

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Target

DAZL

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. IHC: In paraffin embedded Mouse Testis shows cytoplasm staining of spermatogonia. Additional staining is observed in the cytoplasm of elongating spermatids but this additional staining is not in agreement with results obtained by another antibody (Ruggiu et al WB: Approx 35kDa band observed in human ovary and human testes lysates (calculated MW of 33.2kDa according to NP_001342.2). Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence

GNGPQKKSVD

Clonality

Polyclonal

MW

33.2

Entrez

1618

Dilution Range

ELISA: 1:32000, WB: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.


250kDa
150kDa
100kDa

50kDa
37kDa
25kDa
20kDa
15kDa

orb19110 (0.1ug/ml)
staining of human
OV...

250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

Western blot analysis
of Rat Testes lysa...



Immunohistochemical
staining of Mouse
Te...